

The University of Maine

DigitalCommons@UMaine

FAQ & Health Advisory

UMaine COVID-19 Community Archive

1-2020

COVID-19_Provide Saliva for the Shield T3 COVID-19 Test

University of Maine

Follow this and additional works at: https://digitalcommons.library.umaine.edu/c19_health




Part of the [Higher Education Commons](#), [History Commons](#), and the [Medicine and Health Sciences Commons](#)

This Public Health Advisory is brought to you for free and open access by DigitalCommons@UMaine. It has been accepted for inclusion in FAQ & Health Advisory by an authorized administrator of DigitalCommons@UMaine. For more information, please contact um.library.technical.services@maine.edu.



Provide saliva

for the Shield T3 COVID-19 test



1 Remove the cap and hold on to it.




2 Attach the funnel securely to the vial.



3 Deposit 1 ml to 1.5 ml of saliva into the funnel.



4 Remove funnel and screw the cap on the vial. Tighten cap to ensure vial is sealed. Put vial in the plastic bag.



5 Take the bagged vial and funnel to checkout. Hand in vial. Discard funnel in waste container.